

YOU **LECAR**



CATALOGUE2025/26

TABLE OF CONTENTS

YoutecarX base	pg3
Core & System	pg4
Unipolare	pg6
Bipolar	pg8
Neutro Dinamico®	pg9
Automatic Apps	pg10
Joule Mode	pg12
Intimity	pg13
Beauty	pg14
Trigger & Fascia	pg16
Utilities	pg17
Details	pg18
Utilities	pg19



YTX001 - YOUTECAR X BASE

Youtecar X Base, 150W, Unipolar, 500 kHz



CORE&SYSTEM

POWER

200W

200W Power Module

250W

250W Power Module

FREQUENCY

FREQ_STAND

Standard frequency Pack
(500 kHz, 750 kHz, 1MHz)

FREQ_ONDEM

On Demand Frequency Pack
(Standard Frequency Pack
+ one selectable frequency
between 300 kHz and 1 MHz,
in 1 kHz steps)

EMISSION

ATHERMY

HEA (High Energy Athermy),
Hi Energy Pulsed Athermy

ATHERMY-PACK

HEA (High Energy Athermy)
+ SENSITIVE (Soft Emission)

SENSITIVE

SENSITIVE (Soft Emission),
Low Energy Continuous Athermy



CORE&SYSTEM

AUX1

Multifunction AUX Output,
with MULTI Cable included

JOULE-MODE

Additional JOULE working mode,
with measurement of energy delivered



UNIPOLAR

YTX-ELCAP41

Unipolar Capacitive Electrode,
1um gold, diameter 36 mm



YTX-ELCAP61

Unipolar Capacitive Electrode,
1um gold, diameter 56 mm



YTX-ELCAP81

Unipolar Capacitive Electrode,
1um gold, diameter 76 mm



YTX-ELNH

Neutral Electrode Plate in surgical steel
140 x 200 mm



YTX-CAVO-U

RF Flex Cable with UTG Connector
for Unipolar applications



YTX-MAN-U

Handpiece for Unipolar
Resistive and Capacitive applications



UNIPOLAR

YTX-ELRES41

Unipolar Resistive Electrode
in surgical steel, diameter 36 mm



YTX-ELRES61

Unipolar Resistive Electrode
in surgical steel, diameter 56 mm



YTX-ELRES81

Unipolar Resistive Electrode
in surgical steel, diameter 76 mm



BIPOLAR


YTX-ELCAP41B

Bipolar Capacitive Electrode,
1um gold, diameter 36 mm




YTX-ELCAP61B

Bipolar Capacitive Electrode,
1um gold, diameter 56 mm



YTX-ELCAP81B

Bipolar Capacitive Electrode,
1um gold, diameter 76 mm



YTX-LIFT

Dermalift Bipolar Resistive Electrode
with spherical heads in surgical steel



YTX-VERT

Bipolar Paravertebral Resistive Electrode
with heads in surgical steel



YTX-MAN-B

Handpiece for Bipolar Resistive
and Capacitive applications



NEUTRO DINAMICO®



YTX-CAVO-MULTI

RF Multifunction Flex Cable for **Bipolar**, **Neutro Dinamico®** and **Static-Automatic** applications



YTX-MAN-ND

Neutro Dinamico® Handpiece for Unipolar Resistive Electrodes




AUTOMATIC APPS

For application with cable ref YTX-CAVO-MULTI


YTX-CONTOUR/L

Flexible Large Electrode 48x48 mm (4 pcs.)




YTX-CONTOUR/N

Contour Flexible Neutral Electrode 100x81 mm (2 pcs.)



YTX-CONTOUR/S

Contour Flexible Small Electrode 48x48 mm (4 pcs.)



YTX-OCTOPUS

Octopus Kit for Contour Automatic application consisting of:


- 8 Contour Large Kapton Electrodes
- 8 Contour Small Kapton Electrodes
- 2 Contour Neutral Kapton Electrodes
- 1 Hub Box Black
- 1 Hub Box White
- 2 Elastic Bands 100 cm
- 1 Roll of Kinesio Tape
- 12 Conductive Gel adhesive Patches, 40x40 mm
- 12 Conductive Gel adhesive Patches, 50x50 mm

AUTOMATIC APPS

For application with cable ref YTX-CAVO-MULTI



HUB-B

White-Type Contour interconnection Hub



HUB-N

Black-Type Contour interconnection Hub



YTX-ELAUTO

Static-Automatic application Kit consisting of:

- Wiring cable with connections
- 4 Resistive electrodes in surgical steel with clips
- 2 Elastic bands 100 cm



JOULE MODE

JOULE-MODE

Enables the Joule Mode Software

Includes free Octopus Kit for Contour Automatic Application consisting of:

- 8 Contour Large Kapton Electrodes
- 8 Contour Small Kapton Electrodes
- 2 Contour Neutral Kapton Electrodes
- 1 Hub Box Black
- 1 Hub Box White
- 2 Elastic Bands 100 cm
- 1 Roll of Kinesio Tape
- 12 Conductive Gel adhesive Patches, 40x40 mm
- 12 Conductive Gel adhesive Patches, 50x50 mm



INTIMITY

For application with cable ref YTX-CAVO-MULTI

YTX-ENDO

Intimity Handpiece complete with surgical steel
Endocavitary Unipolar Probe, diameter 20 mm



YTX-ENDOA

Intimity Handpiece complete with surgical steel
Endocavitary Unipolar Probe, diameter 12 mm



YTX-ENDOB

Intimity Handpiece complete with surgical steel
Endocavitary Bipolar Probe, diameter 20 mm



YTX-POINT

Unipolar Trigger Points Electrode with spherical tip
in surgical steel



PROBE

Spare Unipolar **Intimity** Probe in surgical steel,
diameter 20 mm (3 pcs.)



PROBE-A

Spare Unipolar **Intimity** Probe in surgical steel,
diameter 12 mm (3 pcs.)



PROBE-B

Spare Bipolar **Intimity** probe Set, in surgical steel and Delrin,
diameter 20 mm (3 pcs.)



BEAUTY

YTX-ELCAP41

Unipolar Capacitive Electrode,
1um gold, diameter 36 mm



YTX-ELRES41

Unipolar Resistive Electrode
in surgical steel, diameter 36 mm



YTX-ELCAP61

Unipolar Capacitive Electrode,
1um gold, diameter 56 mm



YTX-ELRES61

Unipolar Resistive Electrode
in surgical steel, diameter 56 mm



YTX-ELCAP81

Unipolar Capacitive Electrode,
1um gold, diameter 76 mm



YTX-ELRES81

Unipolar Resistive Electrode
in surgical steel, diameter 76 mm



BEAUTY


YTX-ELCAP41B

Bipolar Capacitive Electrode,
1um gold, diameter 36 mm




YTX-ELCAP61B

Bipolar Capacitive Electrode,
1um gold, diameter 56 mm



YTX-ELCAP81B

Bipolar Capacitive Electrode,
1um gold, diameter 76 mm



YTX-MAN-B

Handpiece for Bipolar Resistive
and Capacitive applications



YTX-LIFT

Dermalift Bipolar Resistive Electrode
with spherical heads in surgical steel



YTX-POINT

Unipolar Trigger Points Electrode with
spherical tip in surgical steel




TRIGGER & FASCIA

For application with cable art. YTX-CAVO-MULTI




YTX-IAS TM

Multifunction Fascia Tool with blade in 5 mm thick surgical steel



YTX-POINT

Unipointer Unipolar Trigger Points
Electrode with spherical
tip in surgical steel



UTILITIES



YTX-BTROLLEY

Youtecar X IATA Trolley Bag
with accessories holder



YTX-CARR

Youtecar X complete trolley
with accessories shelf



DETAILS



YTX-BLU

Included with the Basic Version



YTX-RED

YoutecarX with Red finish



YTX-SILVER

YoutecarX with Silver finish



YTX-WHITE

YoutecarX with White finish



YTX-YELLOW

YoutecarX with Yellow finish



UTILITIES

F7820	Neutral Adhesive Gel Electrode 217x105 mm (25 pcs. package)
PG955/10M	Elastic Band for Automatic Electrodes positioning, 10x100 cm
68531710C	Cable with clamp for disposable Neutral Electrode, Black connector
68531710CB	Cable with clamp for disposable Neutral Electrode, White connector
MANH	Neutral Electrode Handpiece in surgical steel, diameter 20 mm
PATCH45	Self-adhesive conductive gel Patch 45 mm (12 pcs. package)
PATCH65	Self-adhesive conductive gel Patch 65 mm (12 pcs. package)
YT-CREAM1000	Yotecar Cream - Conductive Cream 1000 ml (6 pcs. package + 1 dispenser included)
YT-CREAM150	Yotecar Cream - Conductive Cream 150 ml

